



# ONLINE CHALLENGES



VEX VRC HIGH SCHOOL 2022/2023  
Reverse Engineering Online Challenge Report

The Evolution Of Hearing Aids over Time



7282Z

Auto Flashes: Zenith

Franklin Central High School Robotics Club  
6215 S Franklin Rd, Indianapolis, Indiana 46259, United States of America

By Jonathan, Leiam, Armaan, Tory, and Kumuyo

Word count: 436/500

Not including Title page, credits, complete parts lists, image captions, and any appended citations.  
Essentially, all of the "core content" is included in the word count.

<u><b>Table Of Contents</b></u>		<u><b>Page Number</b></u>
1. <i>Introduction</i>		<i>Page 3</i>
1.1 <i>General Timeline</i>		<i>Page 4</i>
1.2: <i>Our Modern Understanding of Hearing</i>		<i>Page 5</i>
1.3 <i>Audiogram Analysis</i>		<i>Page 6</i>
2. <i>Plan of action</i>		<i>Page 7</i>
2.1 <i>Meet the Hearing Aids!</i>		<i>Page 8</i>
3. <i>Hypothesis</i>		<i>Page 8</i>
4. <i>Disassembly</i>		<i>Page 9</i>
4.0.1 <i>Tools Required</i>		<i>Page 9</i>
4.0.2 <i>Custom Vex-IQ Stand</i>		<i>Page 10</i>
4.1 <i>Level 1: Non-Computing Components</i>		<i>Pages 11-12</i>
4.2 <i>Non-electrical Components</i>		<i>Page 13</i>
4.3 <i>Electrical Components</i>		<i>Page 14</i>
5. <i>Parts Lists</i>		<i>Page 15</i>
5.1 <i>Non-Computing Components</i>		<i>Pages 15-21</i>
5.2 <i>Electrical Casing</i>		<i>Pages 22-26</i>
5.3 <i>Electrical Components</i>		<i>Pages 27-35</i>
5.3.1 <i>Attempts to Identify Parts</i>		<i>Page 27</i>
5.3.2 <i>Location of Electrical Components</i>		<i>Page 28</i>
6. <i>Analysis</i>		<i>Page 35</i>
6.1 <i>Weight Distribution</i>		<i>Page 35</i>
6.2 <i>Comparison to Assembled Hearing Aid</i>		<i>Page 36</i>
6.3 <i>Hypothesis Evaluation</i>		<i>Page 37</i>
6.4 <i>Functionality Comparison Across Brands</i>		<i>Page 38</i>
6.5 <i>Power Flow</i>		<i>Page 39</i>
6.6 <i>Theoretical Circuit Schematic of the Moda II</i>		<i>Page 40</i>
6.6.1 <i>Why we Believe the Schematic to be Accurate</i>		<i>Page 41</i>
7. <i>Reassembly Attempts</i>		<i>Page 42</i>
7.1 <i>Electrical → Electrical Casing</i>		<i>Page 42</i>
7.2 <i>Electrical Casing → Complete Hearing Aid</i>		<i>Page 42</i>
8. <i>Conclusions and Lessons Learned</i>		<i>Page 42</i>
9. <i>Works Cited</i>		<i>Page 42</i>
		<i>Pages 43-45</i>

# 1. Introduction

Our Robotics Team Consists of 5 extremely hard workers who all try to make our robot, and team, the best it can be. Despite being in one of the most competitive robotics regions in the world, we've never quit. Instead, we've always searched for how we can improve.



Figure 1.0.1 So far During the 2022-2023 Season, We have won a Design, Innovate, and Create Award!

When we first heard about the Reverse Engineering Challenge, We all had various ideas of what item we should reverse engineer. Out of a Nintendo Wii, First generation Ipod and various other devices, we settled on an older pair of hearing aids. One of our team members, Jonathan Reid, is 40% Deaf and typically wears hearing aids, however no one on our team can confidently explain how they work. As a result, we've decided to discover how these machines that have helped our team member throughout his entire life actually function.




Figure 1.0.2 and 1.0.3: Team Captains Jonathan (left) And Liam (Right) Working on the Robot while Jonathan wears his newest hearing aids.

# 1.1

## How Hearing Aids have Evolved Over Time





Figure 1.1.1: Before starting this challenge, the team wanted to get a grasp on how hearing aids have evolved ever since their conception before delving headfirst into the hearing aids we will disassemble (shown in green on the timeline). We have always believed that understanding the bigger picture of a concept helps with all the details, hence why we've made this timeline. The timeline template is provided by Slidesgo.

The timeline Can be viewed [HERE](#) If one would like to see it in better detail.

# 1.2

Hearing depends on a series of complex steps that change sound waves in the air into electrical signals. Our auditory nerve then carries these signals to the brain. Also available: [Journey of Sound to the Brain](#), an animated video.



Source: NIH/NIDCD

Figure 1.2.1: This Diagram, Provided by the National Institute on Deafness and Other Communication Disorders, shows how the normal ear works. The sound waves travel into the ear where they are processed by the bones and the cochlea into electrical signals for the brain.






Figure 1.2.2: This Graphic Provides a graphic to show how the cochlea gathers sound pitches normally. Although there is no definitive answer, doctors believe that Jonathan's cochleas don't have enough turns in the spiral structure, causing his hearing loss.

# 1.3



Jonathan's Audiogram




Figure 1.3.1: This is a picture of Jonathan's most recent audiogram, which measures his hearing capabilities. The right side of the graphic simplifies the data, showing that Jonathan's hearing falls into the moderate to severe range past 2000 Hertz.

## 2. Plan of Action

Before touching anything, we decided to plan out our process as pictured below




Figure 2.0.1: Process for disassembling and analyzing the hearing aids.

## 2.1

### Jonathan's Hearing Aids Over the Years




Figure 2.2: It should be noted that Phonak And Unitron are sister companies, both owned by Sonova. This means that technology is or can be freely shared between companies, making a viable comparison between the Moda II's and the other two Phonak Models. All manuals, which give vital information about each hearing aid, are linked in the works cited section at the end of this report.

### 3. Hypothesis

If sound waves are around the hearing aid, then the hearing aid will amplify the sound because of the microphone and speaker within it.




Figure 3.1: Our Visual Representation of how we think a hearing aid works. Sound travels through a microphone, and then is amplified and the new sound waves travel through the rest of the system straight into the ear.



## 4. Disassembly

### Pre-disassembly questions

- How does the battery power the System?
- Where/how does sound get amplified or transferred?

#### 4.0.1



Gerber Pliers




Flashlight



"\_\_\_\_\_" Screwdriver



Safety Glasses




Custom Hearing Aid Cleaner



Dial Caliper (for measurements)



15x macro phone lens attachment  
Courtesy of the Franklin Central  
Photography Class



Tweezers

Figure 4.0.1: All of the Tools used To disassemble the hearing aids.

## 4.0.2




Figure 4.0.1: Here is our plan for a custom Stand to display the Hearing aids if/when they are not being disassembled, as only one Modell II hearing Aid will be disassembled. The Parts Library used is free to the Public and is attributed in the sources. You can view the Onshape Model [HERE](#).




Figure 4.0.2 This is How the stand turned out in real life based on the model we created.

# 4.1

## 4.1.1 Moda II Video

Recap:

**How To Disassemble The Non Computing Components of a Phonak Sky-V90 P Hearing Aid**




Figure 4.1.1.1: Make sure that the batteries are out of the system before disassembling anything.

## 4.1.2 Phonak Sky V-90 P Video

Recap:

**How To Dissassemble The Non Computing Components of a Phonak Sky-V90 P Hearing Aid**





Figure 4.1.2.1: The steps of how to disassemble The Phonak Sky V-90 hearing aids.


### 4.1.3 Phonak Naida P Video

Recap:

#### How to Disassemble the Non-computing components of a Phonak Naida P Hearing Aid



Step 1: Pull inner tube out of earmold



Step 2: Pull inner tube and ear hook apart



Step 2: unscrew





Figure 4.1.3.1: How to disassemble the Phonak Naida Hearing aids. There is no battery here, so you just have to verify the hearing aid is off before starting by holding the power button until the indicator light turns red


***From here, we delve into the Moda II hearing aids to see how they operate on a piece-by-piece scale.***

## 4.2


### How to Disassemble The computing Components of a Modal II Hearing Aid




Step 1: Remove Battery Door by pulling the hinge Pin




Step 1: Remove Top Shell with Tweezers.



Step 3: Pry open Sides to loosen the Safety Pin



Step 4: Use Screwdriver to push pin slightly Out



Step 5: Use Pliers to pull pin out safely



Step 6: Push internals out through the back (because the top shell is gone)

Only push it just enough to not damage internal components



Step 7: Use tweezers to pull out the casing the rest of the way




Figure 4.2.1: Our steps for disassembling The hearing aid into the electrical casing.

## 4.3

### How to Disassemble The Electrical Components of A Modall II Hearing Aid

Step 1: Poke Circuitry out of casing

Step 2: There's Pieces in the Bottom, so we must remove those first

Step 3: Use Screw Driver to poke out bottom

Step 4: Pull out speaker and tube connector

Step 5: We also have to remove these Screws

Step 6: Clean gunk off screw

BE CAREFUL!: because these Hearing aids are around 15 years old, The battery door connectors have become rusty, causing it to break

Step 7: Unscrew the screw (both sides)

Step 8: Use screwdriver to create an opening for the microphone

Step 9: **Carefully** Pull out circuitry from Casing

Step 10: Fully take casing off

Step 11: Take off Silicone Grippers on circuitry

Final Circuit


Figure 4.3.1: Our steps for disassembling the electric casing into the circuitry within the hearing aid.

## 5. Parts lists

For Dimensioning, we used the Dial Caliper to measure the farthest distances of the approximate length, height, and depth of each part.

All Measurements are in Inches and labeled by a Length X Height X Depth format

### 5.1

Part Description and Purpose	Image
<p><b><u>Moda II Ear Mold</u></b></p> <p>0.817" X 0.910" X 0.605"</p> <p>Material: Silicone</p> <p>Purpose: Holds the device in your hear/direct sound directly into the ear canal</p> <p>Notable Aspects:</p> <p>Size and exact form are customizable because they Inject your ear with a play-dough type mold to get the size for your ear at the time. This is very similar to die casting for metal, but human friendly.</p> <p><a href="#">Ear molding Putty</a></p> <p><a href="#">Companies like Beltone Create Custom Ear Molds</a></p>	

**Moda II Inner Tubing**

0.915" X 0.550" X 0.175"

Material: Silicone

Purpose: Direct Sound and connect earmold with Ear hook

Notable aspects: Much more flexible which allows for greater comfort but also an easier chance of tearing.



**Moda II Ear Hook**

0.885" X 0.310" X 0.210"

Material: Silicone

Purpose: Connect hearing aid and tubing and provide more structure along the curvature of the ear

Notable Aspects: Has a much more sizeable attachment to the hearing aid unit than the Phonak brand ones





**(Assembled) Moda II hearing aid unit**

Battery Door Closed

1.125" X 0.494" X 0.260"

Battery Door open

1.543" X 0.650" X 0.260"

Material:

It's... complicated (because it has so many parts within it, which you will see later)



Purpose: Receives and amplifies sound.

Notable Aspects

Product has been discontinued, but can be found for resale on places like [ebay](http://ebay.com).

Data Sheet:

<https://www.unitron.com/content/dam/unitron-2014/documents/english/Older%20Products/Next/specsheet/english-next-specsheet-letter-ne-modaii.pdf>



**Size 312 Battery**

0.145" X 0.305" diameter

Material: Zinc

Purpose: Powers the system

Notable Aspects:

[Here's how batteries work](#)  
[More Information](#)

Much smaller than the size 13 battery.  
(It is approximately 3/4th it's height in fact)

[Available online](#)



**Phonak Sky v90 P Ear Mold**


0.877" X 0.135" X 0.4"

Material: Silicone

Purpose:  
Hold the device in your hear/direct sound directly into the ear canal

Notable Aspects:

Clear material allows the earmold to blend in more with one's natural ear so that they don't stand out as much.




**Phonak Sky v90 P Inner Tube**

0.815" X 0.421" X 0.175"

Material: Silicone

Purpose: Direct Sound and connect earmold with Ear hook.

Notable Aspects: Much firmer structure than either of the 2 other inner tubes... We think this is because of a thicker silicone or even flexible plastic layer.



**Phonak Sky v90 P Ear Hook**

0.820" X 0.345" X 0.175"

Material: Hard Plastic of some sort

Purpose:

Connect hearing aid and tubing and provide more structure along the curvature of the ear.

Notable Aspects: Has threading on it for tighter grip with the hearing aid unit.

[Available Online](#)



**Phonak Sky v90 P Hearing Aid Unit**

Door Closed  
1.31 l, 0.581 h, 0.335 d

Door Opened  
1.742" X 0.762" X 0.335"

Material: Various

Purpose: Receives and amplifies sound

Notable Aspects:

Has a volume control switch on the top, as well as very durable casing

[Website overview of a newer innovation of this model, The Sky Marvels](#)

Data Sheet:  
[https://www.phonakpro.com/content/dam/phonakpro/gc\\_hq/en/products\\_solutions/hearing\\_aid/sky\\_v/documents/datash eet\\_sky\\_v\\_sp.pdf](https://www.phonakpro.com/content/dam/phonakpro/gc_hq/en/products_solutions/hearing_aid/sky_v/documents/datash eet_sky_v_sp.pdf)



**Size P13 Battery**

0.205" X 0.305" Diameter

Material: Zinc

Purpose: Powers the system

Notable Aspects:

[Here's how batteries work](#)  
[More Information](#)

Larger than size 312 batteries by about 0.05 inches, but comparatively it makes a big difference. This battery is about 25% taller than the size 312 Battery used with the Moda II hearing aid.

[Available Online](#)



**Phonak Naida Earmold**


0.755" X 0.943" X 0.591"

Material: Silicone

Purpose:  
Hold the device in your hear/direct sound directly into the ear canal

Notable Aspects:

This is a skeleton Earmold, which has a hole in it that allows for more breathability, but is also more fragile.




**Phonak Naida Inner Tubing**

0.841" X 0.54" X 0.200"

Material: Silicone and metal (of some sort) for the bearing

Purpose: Direct Sound and connect earmold with Ear hook

Notable Aspects: Instead of a plastic fitting, this inner tube uses a metal bearing to firmly fit inside of the earmold



**Phonak Naida Ear Hook**


0.800" X 0.400" X 0.210"

Material: Hard Plastic of some sort

Purpose: Connect hearing aid and tubing and provide more structure along the curvature of the ear

Notable Aspects:

Has both threaded attachment and a bearing type fitting



**Phonak Naida Hearing Aid Unit**

1.415” X 1.142” X 0.335”

Material: Various

Purpose: Receives and amplifies sound

Notable Aspects:

No battery, as it is rechargeable. It also has Bluetooth support as well as its own app you can set custom settings.

[Website Overview](#)

Data Sheet:

<https://drive.google.com/file/d/1UngowFD9Bq56uzjGKNr6Cs0Eiyz8soeB/view?usp=sharing>



## 5.2









Figure 5.2.1 and 5.2.2: This is a complete render and an exploded drawing view of the digital recreation of the Moda II hearing aid using Onshape. While it's not as intricate as the real hearing aid, this model provides an accurate representation and visualization on how the parts fit together. [VIEW IT HERE](#)

Part Description	Image	Recreated Model Drawing Of Part
<p><b><u>Brown Cover Shell</u></b></p> <p>1.004" X 0.455" X 0.250"</p> <p>Material: Hard Plastic</p> <p>Purpose:</p> <p>Hold white casing securely.</p> <p>Notable Aspects:</p> <p>Might have been made using injection molding, based upon a seam we noticed inside the casing.</p>		
<p><b><u>Top Slot Protector</u></b></p> <p>0.515" X 0.200" X 0.106"</p> <p>Material: hard plastic</p> <p>Plastic</p> <p>*we believe it to be polystyrene, or plastic type 6 based on the number 6 present on the bottom along with material testing*</p> <p><a href="#">Plastic types</a></p> <p>Purpose:</p> <p>Protects the microphone and provides a filter so debris doesn't get in</p> <p>Notable Aspects:</p> <p>Made a different color than brown casing to emphasize it</p>		

**White Casing Unit**

0.71" X 0.585" X 0.175"



Material: Plastic

Purpose: Hold Circuitry

Notable Aspects:

Serial number of the hearing aid is visible here. However, the hearing aids have been discontinued since sometime in the 2000s. That, combined with patent laws, makes it near impossible to systematically identify the parts in terms of the official company unless they imported it from someone else. We even reached out to Unitron, the manufacturer, to no avail.

The Serial number is 0914XOTE8



**Battery Door**


0.645" X 0.388" X 0.246"

Material: hard plastic  
\*potentially type-6 but no confirmation\*

Purpose:  
Holds battery and compressed it when closed that powers the system.

Notable Aspects:

The top notches at the top plug into the battery plug which helps power the system.





**Hinge Pins**


0.236" X 0.027" diameter

Material: metal

Purpose:  
Hold battery door and white casing

Notable Aspects:

Tolerance in the hole and pin allows pieces to swivel while being held in tight by the casing.



**Bottom Silicone Holder**


0.36" X 0.025" X .17" X 0.165"

Material: Silicone

Purpose:  
Protect Speaker and Pipe Plug.

Notable Aspects:

Held onto the white casing simply by tape.



**Silicone Circuit Holders**

Right Microphone/processor

0.145" X 0.185" X 0.1"

Left Microphone/processor holder




0.143" X 0.145" X 0.125"

\*there was another part identical to the one on the right, but it has been misplaced, however, it is evident on the model\*

Material: Silicone

Purpose: Holds the microphone/amplifier in a raised position in the white casing

Notable Aspects:  
Incredibly Fragile and susceptible to tearing.



## 5.3

### 5.3.1

#### Reverse Engineering Project Help/Inquiry >



Jonathan Reid

to mediarelations, mediarelations

Thu, Dec 22, 2022, 11:33 AM (13 days ago)



Hello, My name is Jonathan Reid and I am a Senior at Franklin Central High School in Indianapolis Indiana. More importantly, I am a head leader of the Vex Robotics team offered at my school. As a result of being in this program, members are eligible to participate in a variety of "challenges" for a chance to qualify for the International or Worlds event held in May of 2023.

One of these challenges is the [vrc-high-school-reverse-engineering-challenge](#), which asks you to reverse engineer an electronic device to discover how it works. Because of my lifelong hearing loss, I've chosen to reverse engineer my first pair of hearing aids, The [Next Essential Moda II BTE](#), as manufactured by Unitron.

However, despite my best attempts at researching the individual parts found, I simply couldn't find an answer to everything. The only electrical piece I was able to find online was [This speaker unit used](#) in the hearing aid.

***The Current Draft Is Linked or Attached at the Bottom of This Message***

If there is any way Unitron or Sonova as a whole could help me in identifying these parts I would be forever thankful, but I understand if it's not feasible because of various Patents or Legal Limitations.

If there's something, **anything**, that can be done to help, please reply to this email as soon as possible.



***We tried directly reaching out to Unitron and using a microscope to find part numbers or identifiers for parts of the hearing aid. Both attempts yielded little success due to either company trademarks or the fact that these circuits are custom made for this discontinued brand of hearing aid.***

### 5.3.2

## Location of major electrical components of Moda II hearing Aid






Figure 5.3.2: The locations of all visible components within the Hearing Aid Casing to get a sense a relativity and how each piece works off the others

### 5.3.3



Snippet of our electrical section of our CAD recreation of the electrical components found [HERE](#)



Number	Part
1	Ribbon Cables
2	Amplifiers
3	Sound Transfer Pipe
4	Battery Plug-In
5	Microprocessor
6	Battery Holder
7	Screw
8	Speaker
9	Pipe Plug
(not pictured)	Wiring

Figure 5.3.2.1 and Figure 5.3.2.2: Like Before with the electrical casing, here is our rendered recreation of the electrical components as well as the Exploded view and Bill of Materials. [VIEW IT HERE](#).

**Part Description**

**Image**

**Recreated Model of Part**

**Ribbon Cables**

0.005" X 0.136"

Length is unable to accurately determine because of odd shaping. Approximated length is about 0.84 when attempted to lay flat

Material:



[Ribbon Cables and what They are made of](#)

Purpose:

Transfers power from battery system to amplifier system

Notable Aspects:

Appears to have three copper strands per cable



**Sound Processors/Amplifiers**

0.288" X 0.089" X 0.157"

Material: Metallic

Purpose:

Takes power and likely transfers and amplifies the pitches that the user cannot hear


Notable Aspects:

We found two inscriptions on the amplifiers, but everywhere we have looked has led to a dead end even after months of research. The inscriptions read [5208MB](#) and [26403CX](#) respectively.

[Amplifier Explanation 1](#)

[Amplifier Explanation 2](#)

Here are some examples of amplifiers we could find online: [1](#) [2](#)



**Microphone/ Sound Transfer Pipe**

0.245" X 0.045" diameter


Material: Metallic

(we want to say copper because of its color and texture, but cannot confirm)

Purpose:

Transfer base pitches back into the system that doesn't need processing from the amplifiers.

The wires were attached by some bonding material, as marks are left where they were



**Battery Plug In**

0.109" X 0.110"



Material: plastic casing with metal prongs.

Purpose:

Receive power from battery door.

Notable Aspects:

We originally thought this to be the microphone, but after further inspection with the completed right hearing aid, it clearly connects to the battery door and adds power to the system.



**Microprocessor**

0.210" X 0.070" X 0.151"

Material:

[What a microprocessor is made out of](#)

Purpose:



Performs logic (I.E. what pitches are heard and which ones need to be adjusted)

Notable Aspects:


[Further Information about microprocessors](#)

[Various Texas Instrument Amplifiers Microprocessors](#)

[Shopping link for Microprocessors](#)



*This image was taken through a microscope for a better look at the circuitry that is too small to dissect.*



**Battery Connectors (Metal Plate)**

0.241" X 0.088" X 0.005"

Material:


Metal  
Aluminum, Iron, or steel perhaps because of the conductive properties

Purpose:

Hold battery and help direct charge

Notable Aspects:

Fragile because of rusting from age





**Screws**


0.092" Total Length  
0.046" Screw Head Diameter  
0.027" Thread Diameter

Material:  
Metallic (Aluminum, Iron, or steel)

Purpose:  
Hold Battery Connectors in tight

Notable Aspects:

Susceptible to Rust



**Speaker Unit: Sonion Speaker**  
26J02/7

0.285" X 0.113" X 0.13"

Material(s): Metal and circuitry




Purpose:  
Amplify the volume of the sound and then send it to the "pipe plug"

Notable Aspects:

This part was independently bought by Unitron from Sonion for their product so we were able to find it online.

DATA SHEET:

[https://www.sonion.com/wp-content/uploads/26J027-2003827.pdf?\\_gl=1\\*\\_rvzh3d\\*\\_ga\\*MTA3MDQ1OTEzLjE2NzEzOTQ1MDM.\\*\\_up\\*MQ..](https://www.sonion.com/wp-content/uploads/26J027-2003827.pdf?_gl=1*_rvzh3d*_ga*MTA3MDQ1OTEzLjE2NzEzOTQ1MDM.*_up*MQ..)




**Pipe Plug**

0.267" in Length  
0.138" ring diameter  
0.070" metal diameter  
0.067" silicone diameter

Material: silicone and metal

Purpose:

Takes Speaker Noise and transfers it to Hear hook



**Wiring**

0.002" d

Can't find the length of wire using dial caliper

Material: small refined copper or any other conductive metal


[what wires are made out of](#)

Purpose:

Transfer sound and electricity

Notable Aspects:

Labeled by a red wire and a green wire, but we can't tell why there's a distinction. However, we believe it may be because one wire transfers the "base" sound pitches while the other transfers the "altered" pitches



Model Not Available because of sizing and flexibility issues with extrusions

## 6. Analysis

### 6.1

# Weight Distribution Chart (in Grams)




Figure 6.1.1: This a Pie chart showing the weight distributions of each hearing aid as a whole, marked by the color. However,we subsided each chart with all of the disassembled components we found. For example, the Earmold for the Phonak Naida hearing aid was about 12% of the total mass of all of the hearing aids, but only about a 4th of the weight of the Phonak Naida itself.

## 6.2

Key:

Battery Connector

Battery Plug

Amplifiers (visible through small vent)

Pipe Plug

Speaker

Motherboard









Figure 6.3.1: This is a diagram showing how one can find the electrical pieces we found in our disassembly in an assembled version of the hearing aid. You can clearly see the battery connector, battery plug, and pipe plugs without touching anything. We can also see/infer the location of other pieces, such as the amplifiers, through vent holes protected by the top slot protector

# 6.3

Hypothesis: If sound is around the hearing aid, then the hearing aid will amplify the sound because of the microphone and speaker within it.

Actual: If sound waves are picked up by the microphone, then the hearing aid will process and amplify certain un-hearable pitches (as defined by the audiologist for the patient) and then send all sounds received out directly into the users ear because of the Tubing systems



## What we got Right

## What we Got Wrong

There is Indeed a Microphone, but just not where we thought it was	The device on top that we thought was a microphone is actually where the battery door connects to the system
The Speaker System is on the bottom of the hearing aid and does send sound through the earhook and tubing	Sound Pitches themselves are actually amplified because of the Processors in (in relation to the hypothesis picture, the top right of the hearing Aid)
General Pathing of how Sound was Correct, albeit vastly oversimplified	The battery plays a much larger role than expected as shown through its multitude of attachments to the circuit
	The circuit itself is far more complicated than hypothesized

Brown Case is Translucent for detailing

Figure 6.4.1: Here is an analysis we created to contrast our hypothesis with our actual findings throughout this reverse engineering project.

## 6.4




Category	Moda II	Phonak Sky	Phonak Naida
Picture			
Overall Sizing (inches)	2" X 1 ¼" X 1 ½"	1 ½" X 2" X ¾ "	2 ½ " X 2" X ¾"
Weight (in Grams)	4.2	5.4	5.9
Power Sources			
Average lifespan (before replacing and or recharging)	3- 10 days	6- 14 days.	24 hours of continuous use in a single charge
<u>Battery?</u>			
312			
<u>P 13</u>			
Rechargeable?			
Charge Time	Instant	Instant	
Non Computing Parts			
<u>Earmold</u>			
Skeleton Structure?			
Clear?			
<u>Ear Hook</u>			
Threaded Attachment To hearing Aid?			
The *Actual* Hearing Aids Themselves			
Volume Switch?			
Bluetooth connectivity?			
Water Proofing?			

Figure 6.4.1: This is our visual representation of the comparisons we could make between all of the hearing aids Jonathan has owned. The purple cells are ones with variable data for each hearing aid. The blue cells are used as titles and for signifying a hearing aid doesn't have a function, where as the green cells signify that they do.

## 6.5

### Power Flow Chart of Moda II Hearing Aid




Figure 6.5.1: This is how we believe this hearing aid works with the newfound knowledge of all the components inside of it.

## 6.6

### What We Think is Happening in the Unitron Next E Moda II Hearing aids *Electrically*





Figure 6.6: this is what we believe a vastly simplified version of the hearing aid would look like on an electrical schematic based on all of our research. The battery would feed directly into the microprocessor, battery plug, and microphone, then split off into various parallel circuits for each of the amplifying tracks. After this, all the connected sound would pass out of the speaker and the remaining electrical signals would go back into the battery, completing the circuit.




## 6.6.1


### Why We Believe This is What the Circuitry looks like



Hearing Aid Schematic with 2X 1.5 volt batteries



An "inexpensive hearing aid circuit that uses just four transistors and a few passive components."



Hearing Aid using audio amplifier LM386



Hearing aid recreated on a breadboard by NotionSaturday

Fig 6.6.1: This is the culmination of the research we used as both an inspiration and way to check our work for our own custom hearing aid circuit schematic. Many of these are made to be vastly simplified versions of hearing aids, but some, such as the prototype made on the bottom right by [NotionSunday](#) shows just how simple a hearing aid could be boiled down too using a breadboard and basic electronic materials.

## 7. Reassembly Attempt

### 7.1 Electrical--> Electrical Casing

Video 7.1: This is our attempt at reassembling the electrical aspects of the hearing aid into the white casing unit. We were largely successful however because of damage from the disassembly, age, or even just tape or glue, it was very hard and potentially dangerous to shove the electricals back into the brown casing...

### 7.2 Electrical Casing ---> Complete Hearing Aid

Video 7.2: ... so we created a second video with the other identical hearing aid to show how that process would ideally work as the electronics were not tampered with on that side.

## 8. Conclusion

Through this project, we've learned how to safely disassemble, analyze, compare, and reassemble tiny electronic devices. On top of that, we've uncovered a bit of biology and history by analyzing how hearing aids as a whole have evolved over time. We also deepened our understanding of concepts such as CAD and the Engineering design process as a whole by trying to understand how these hearing aids were designed by Unitron. Most importantly, we learned how technology provides invaluable support for those put into a natural disadvantage, helping them to live their lives to the fullest.



Figure 8.1: All the colors offered for Jonathan's newest hearing aids, the Phonak Naida Paradise

### What we've learned:

- Even seemingly simple machines are very intricate
- Small products are often the most compact
- It's crucial to document all steps of a disassembly to track pieces and have an easier reassembly process
- Dial Calipers, while precise, do not measure curvatures nor complex shapes accurately
- Hearing Aids have been innovated upon for centuries
- What A Cochlea is and how it works
- Companies can be stingy with their product to protect trademarks
- How to write and read circuit schematics

## Special Attributions/Works Cited

### User Manuals

- Yumpu.com. (n.d.). *Unitron hearing - Moda II - BTE user guide - English - PDF*. Yumpu.Com.  
Retrieved December 29, 2022, from  
<https://www.yumpu.com/en/document/read/19367948/unitron-hearing-moda-ii-bte-user-guide-english-pdf>
- AG, P. (n.d.-a). *User guide sky V*.  
[https://www.phonak.com/content/dam/phonak/HQ/en/solution/products/sky\\_v/documents/user\\_guide\\_sky\\_v\\_029-0370.pdf](https://www.phonak.com/content/dam/phonak/HQ/en/solution/products/sky_v/documents/user_guide_sky_v_029-0370.pdf)
- Phonak. (n.d.-b). *User Guide Naida P-PR (pdf)*  
[https://www.phonak.com/content/dam/celum/User\\_guide\\_Naida\\_P-PR/PH\\_User\\_guide\\_Naida\\_P-PR\\_92x125\\_EN\\_V1.00\\_029-0872-02.pdf](https://www.phonak.com/content/dam/celum/User_guide_Naida_P-PR/PH_User_guide_Naida_P-PR_92x125_EN_V1.00_029-0872-02.pdf)

### Whisper Hearing Centers for Proctoring the Audiogram Featured in 1.2

- Gurnick, S., Au. D. (2019, November 14). *Whisper Hearing Centers*. Whisper Hearing Centers.  
<https://whisperhearingcenters.com/>

### Other References

- Acousticon model A-120 “constellation” vacuum tube hearing aid*. (n.d.). Retrieved December 29, 2022, from  
<https://hearinglosshelp.com/museum/acousticon-model-a-120-constellation-vacuum-tube-hearing-aid/acousticon120/>
- admin. (2020, October 19). *Who invented the hearing aid?* Beck Hearing Aid Centre.  
<https://beckhearingaids.com/who-invented-the-hearing-aid/>

Amg. (n.d.). Electrical & Electronic Symbols. Amplifier Circuit Symbols. Retrieved January 6, 2023, from <https://www.electrical-symbols.com/electric-electronic-symbols/amplifier-circuits-symbols.htm>

Bhatt, A. (n.d.). Hearing aid using audio amplifier LM386. Engineers Garage. Retrieved January 6, 2023, from <https://www.engineersgarage.com/hearing-aid-using-audio-amplifier-lm386/>

*Bluetooth hearing aids are ready for prime time.* (2020, February 25). Bluetooth® Technology Website. <https://www.bluetooth.com/blog/bluetooth-hearing-aids-are-ready-for-prime-time/digital-hearing-aids-New-Generation-Hearing-Blog>

*digital hearing aids - New Generation Hearing Blog.* (n.d.). New Generation Hearing Blog. Retrieved December 29, 2022, from <https://newgenerationhearing.wordpress.com/tag/digital-hearing-aids/>

*ELI: Today's standard for bluetooth hearing instrument communication.* (n.d.). AudiologyOnline. Retrieved December 29, 2022, from <https://www.audiologyonline.com/releases/eli-today-s-standard-for-5102>

Foundation, H. H. (2021, March 4). Hearing aid history: Ear trumpets, European royalty, & earbuds – Hearing Health Foundation. *Hearing Health Foundation*. <https://hearinghealthfoundation.org/blogs/hearing-aid-history-ear-trumpets-european-royalty-earbuds>

GaryC. (2022, June 20). Hearing Aids Circuit. ElectroSchematics.com. Retrieved January 6, 2023, from <https://www.electroschematics.com/digital-hearing-aids-circuit/>

Hearing aid batteries. Alapa'i Audiology - Honolulu Hearing Aids and Hearing Services. (n.d.). Retrieved January 6, 2023, from <https://www.alapaihearingaids.com/hearing-aid-batteries>

Hearing aid. Hearing Aid circuit diagram and instructions. (n.d.). Retrieved January 6, 2023, from <http://www.hobby-circuits.com/circuits/audio/audio-amplifier/759/hearing-aid>

*Home.* (n.d.). Unitron.Com. Retrieved December 29, 2022, from [https://www.unitron.com/us/en\\_us.html](https://www.unitron.com/us/en_us.html)

*How do we hear?* (n.d.). NIDCD. Retrieved December 29, 2022, from <https://www.nidcd.nih.gov/health/how-do-we-hear>

*How to read an audiogram and determine degrees of hearing loss.* (n.d.). Retrieved December

29, 2022, from <https://www.nationalhearingtest.org/wordpress/?p=786>

Human ear. (1998, September 16). *Encyclopedia Britannica*.  
<http://britannica.com/science/ear/Transmission-of-sound-within-the-inner-ear>

*Introducing Phonak Sky V - The new pediatrics product portfolio*. (n.d.). AudiologyOnline.  
Retrieved December 29, 2022, from  
<https://www.audiologyonline.com/interviews/introducing-phonak-sky-v-226-16883>

Iris, R. (n.d.). *Phonak media release: Naida Marvel*.

Josephson, M. (1998, July 20). Thomas Edison. *Encyclopedia Britannica*.  
<http://britannica.com/biography/Thomas-Edison>

Seelisch, A. (2016, August 23). *Hearing solutions*. Hearing Aid Service - Hearingsolutions.Ca.  
<https://www.hearingsolutions.ca/the-history-and-evolution-of-the-hearing-aid/>

Seniorliving.org. (2019, July 25). *Phonak vs Unitron hearing aids - Which hearing aid should you buy?* SeniorLiving.Org. <https://www.seniorliving.org/hearing-aids/phonak-vs-unitron/>

*Sitemap*. (2021, March 25). Resources on Disability Assistance: Your Rights and Benefits.  
[http://alphaonenow.com/userfiles/image/aud\\_audiogram.gif](http://alphaonenow.com/userfiles/image/aud_audiogram.gif)

Sonova. (n.d.). Phonak. Retrieved December 29, 2022, from  
<https://www.phonak.com/en-us?null>

*Starkey Laboratories debuts bluetooth® device for the hearing impaired*. (n.d.).  
AudiologyOnline. Retrieved December 29, 2022, from  
<https://www.audiologyonline.com/releases/starkey-laboratories-debuts-bluetooth-4978>

YouTube. (2016, October 14). DIY low cost stereo hearing aid | hearing amplifier. YouTube.  
Retrieved January 6, 2023, from <https://www.youtube.com/watch?v=bAZbHFOIRfg>